

HEIGHT ADJUSTABLE
180mm



**NEW
PRODUCT
RELEASE!**

**REMOVABLE
ARMRESTS / WINGS**
For SIDE TRANSFERS

**REMOVABLE
LITHIUM-ION
BATTERY**
No Cables

MOBILE
Ø75mm Castors



IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.
The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

ULTRACARE MOBILE LIFT CHAIR

LC0901

FEATURES

- Lift / Rest / Recline armchair with height adjustment.
- Mobile lift chair with removable Lithium-ion battery.
- Ø75mm castors with push handle at rear.
- Handset controls with large buttons.
- Removable armrests and wings for side transfers.
- Incontinent PU backrest / seat / legrest coverings
- Vinyl side coverings with storage pockets either side.
- Added comfort of a pocket spring seat base.
- Quality German actuators.
- Audible and visual warning when battery requires recharging.
- Warranty: Frame: 5 years / Electrics: 2 years / Fabric: 2 years

SPECIFICATIONS All dimensions shown are nominal

Overall: **760mm (W) x 850mm (D) x 1140mm (H)**
Seat: **540mm (W) x 510mm (D) x 530 - 710mm (H)**
Armrest: **100mm (W) x 450mm (L) x 180mm (H)**
Backrest: **540mm int. (W) x 700mm (H): 110 - 145°** (Relevant to seat)
Legrest: **540mm (W) x 450mm (L): 90 - 15°** (Relevant to seat)
Castors: **Ø75mm - Swivel front / Braked rear**
Distance from wall: **600mm**
Product Weight: **65kg**
Maximum User Weight: **160kg**
Power Source: **Removable rechargeable Lithium-ion battery pack**
25.2V - 2550mA - 65Wh

Please contact **REDGUM Brand**
if you require further information.

P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand - 1 Business way, Malaga, Western Australia 6090

redgumbrand.com.au

Technical Specifications

Surface Composition:	100% Polyurethane	
Base Composition:	100% Polyester	
Tensile Strength: EN ISO 1421:	Warp 568 N 72.8% 620 N	Weft 103.2%
Tear Growth Resistance: EN ISO 4674 - 1	Warp 63.8N	Weft 61.9N
Permanent Creasing Behavior (Bally Test): DIN EN ISO 32100/DIN 53359 - 100.000 / cycles	PASS	
Colour Fastness to Rubbing: EN ISO 105-X12	Dry 5	Wet 5
Colour Fastness to Artificial Light: EN ISO 105-B02	>=6	
Resistance to Urine (Urea Solution):	PASS	
Resistance to Hospital-Grade Disinfectant Cleaner:	PASS	

Application & Features



High Scratch Resistance



Resistance to Hydrolysis (10 years)



Durability Performance/ Abrasion Test: > 500.000 cycles (Heavy Use)



Fire proof norms

BS5852 / BS 5852 Part 1 - Cigarette & Match Tests - PASS

EN1021 / EN1021-1/2 - Cigarette & Match Tests - PASS

TB - 117 / CA Tech. Bull. 117 / NFPA 260 - Cigarette Test - PASS



Chemical Solvents-Free



Respects the environment

Clean and Care recommendations



WASH



BLEACH



DRY



DRY



DRYCLEAN



IRON

- Frequent wash may affect the performance of the article.
- To remove dust use a vacuum cleaner.
- Any stain, spills or soiling should be cleaned up promptly with a neutral soap solution using a cloth in smooth circular movements.
- We do not accept responsibility for any damage resulting from inappropriate use or chemical cleaning.