

# MOVO FLIP UP CANOPY COMPATIBILITY CHART

## FOR THE SHOPRIDER RANGE

### IMPORTANT NOTE:

Mobility Scooter must have the following in order to fit a Movo Canopy:

- Socket under the seat that can take a 25 x 25mm square tube
- Overall width of armrests must be less than 630mm
- Distance from bottom of socket to top of backrest (not including a headrest) must be more than 600mm

							
<b>MODEL</b>	<b>ROCKY 8</b>	<b>ROCKY 4</b>	<b>SEKA</b>	<b>889SL</b>	<b>ALLROUNDER</b>	<b>EXPLORER</b>	<b>VENTURER</b>
<b>PRODUCT CODE</b>	889AE	889XL	889ASN	889SL	9LSP	888SE	888IX
<b>MOVO COMPATIBILITY</b>	<b>YES</b>	<b>NO</b>	<b>YES</b>	<b>NO</b>	<b>NO</b>	<b>YES</b>	<b>YES</b>
							
<b>MODEL</b>	<b>ROCKY GOLF</b>	<b>889D</b>	<b>ROCKY 12</b>	<b>VIKING</b>	<b>778SL</b>	<b>778HD</b>	<b>CROSSOVER</b>
<b>PRODUCT CODE</b>	889RGB	889D	889WB	779XLS	778SL	778HD	GK10
<b>MOVO COMPATIBILITY</b>	<b>YES</b>	<b>NO</b>	<b>NO</b>	<b>NO</b>	<b>NO</b>	<b>NO</b>	<b>YES</b>